



# RETIKARE I/P

AIR INJECTOR / GAS FORCED INFUSION

The RETIKARE I/P unit allows you to use as Air Injector Function or Infusion Pressure (Air Forced Infusion) with Temponade for both Anterior/Posterior Surgical Procedures

RETIKARE I/P provides you the following modes of operations.

- A. Air Injector Function (or)
- B. Air Forced Infusion Pressure with Temponade Function

## FEATURES & SPECIFICATIONS

- Positive Air Pressure Output upto 100mmHg
- Inbuild Air Pump (Eliminates External Air Pressure source)
- Temponade Function enable/disable through separate Foot switch
- Used for both Anterior/Posterior Surgical Procedures
- Easy to read 7 segment display
- Feather touch control keys for cut rate Up/Down
- Light weight and Portable
- Electrical power: AC 100-240 Volts, 50/60 Hz, Single Phase
- Weight: 1.2 Kg
- Dimensions: L 225mm x D 199 mm X H 72 mm



Manufactured by:



ICON MEDICARE  
9, HIG-II, TNHB Phase VI, Rainbow Garden,  
Bagalur Road, Hosur - 635 109. Tamil Nadu. INDIA  
PH : +91 4344 329890  
Email: [info@iconmedicare.com](mailto:info@iconmedicare.com) [www.iconmedicare.com](http://www.iconmedicare.com)